

Sales Information

Key Words:

HILIC YMC
Polar Analytes
Scale-Up

No. 110-5 e

YMC
EUROPE GMBH

HILIC with YMC



To allow you to succeed in the analysis of extremely polar compounds, YMC has developed several selectivities, which enables reversed phase work with aqueous eluents: **YMC-Pack ODS-AQ** and **YMC Hydrosphere C18**, (the latter is based on YMC's ultra pure **ProFamily** silica and therefore allows excellent peak shape to be obtained.) If the retention of the polar compounds is insufficient on these C18 phases, **Hydrophilic Interaction Chromatography HILIC** can provide an alternative separation mechanism. YMC offers a variety of very rugged selectivities, which complement the above reversed-phase materials.

• **Silica** • **PVA-Sil** • **Amino** • **Polyamine II** • **Cyano** • **Diol**

The mainly normal phase materials cover a broad range of polarity so that even extremely polar compounds can be retained to provide an optimal analysis. In addition, **HILIC** conditions of high organic phase concentration mobile phases are very compatible with LC/MS systems.



With the wide variety of particle sizes (from 3 µm to 50 µm) for YMC phases, scale-up to prep and process scale is simplified. In addition, the YMC bulk materials can be produced in multi-ton quantities to allow full scale production processes to be established.

Applications:



1. Nicotinic acid
2. Erythorbic acid
3. L-Ascorbic acid

Column: YMC-Pack NH₂ (5 µm, 12 nm)
250 x 4.6 mm ID
Eluent: Acetonitrile / 50 mM NH₄H₂PO₄ (70:30)
Flow rate: 1.0 ml/min
Temperature: 40 °C
Detection: UV at 250 nm, 0.16 AUFS
Injection: 10 µl (0.05-0.1 mg/ml)



1. Panthenol
2. Nicotinic acid
3. Urea
4. Allantoin
5. L-Ascorbic acid
6. Pyridoxine HCl (Vitamin B₆)

Column: YMC-Pack Polyamine II (5 µm, 12 nm)
250 x 4.6 mm ID
Eluent: 50 mM NH₄H₂PO₄-H₃PO₄ (pH = 2.8) / Acetonitrile (20:80)
Flow rate: 1.0 ml/min
Temperature: 37 °C
Detection: UV at 210 nm, 0.16 AUFS
Injection: 10 µl (0.01-5.0 mg/ml)

Sales Information


Key Words:

**HILIC YMC
Polar Analytes
Scale-Up**


No. 110-5 e

YMC
EUROPE GMBH

Applications:



Column: PVA-Sil 100 x 3.0 mm ID
Eluent: Isocratic 95% A / 5% B
A: Acetonitrile / B: Ammonium acetate
200 mM, pH 5,5
Flow rate: 0.9 ml/min



Column: Polyamine II 125 x 4.0 mm ID
Eluent: Gradient 100% A -> 20 min -> 80% A
A: Acetonitrile / B: Ammonium acetate
200 mM, pH 5,5
Flow rate: 0.9 ml/min

Ordering Information

Phase	Pore Size	Particle Size	Order-No.
PVA-Sil	12 nm	5 µm	PV12S05...
Polyamine II	12 nm	5 µm	PB12S05...
Amino	12 nm	3 µm	NH12S03...
Amino	12 nm	5 µm	NH12S05...
Cyano	12 nm	3 µm	CN12S03...
Cyano	12 nm	5 µm	CN12S05...
Cyano	30 nm	5 µm	CN30S05...
Diol	6 nm	5 µm	DL06S05...
Diol	12 nm	5 µm	DL12S05...
Diol	20 nm	5 µm	DL20S05...
Diol	30 nm	5 µm	DL30S05...

Phase	Pore Size	Particle Size	Order-No.
Silica	6 nm	3 µm	SL06S03...
Silica	6 nm	5 µm	SL06S05...
Silica	12 nm	3 µm	SL12S03...
Silica	12 nm	5 µm	SL12S05...
Silica	20 nm	3 µm	SL20S03...
Silica	20 nm	5 µm	SL20S05...
Silica	30 nm	5 µm	SL30S05...

**Simply add the column dimensions to the order numbers, with the length in cm and the column i.d. in two digit mm e. g.:
4.0 mm = 04**